

REF 510-2	6x3mL	EXP 2026-07-22	LOT 30A891C	更新日期： 2021-09 1 / 2
------------------	-------	-----------------------	--------------------	------------------------

批号：30A891C				
项目\仪器\试剂方法	单位	均值	+ / - 2 SD	
(AFP)甲胎蛋白				
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	72.8	55.3-90.3	
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	76.1	57.8-94.4	
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	56.9	43.2-70.6	
(CA 125)糖类抗原 125				
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	268	204-332	
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	241	183-299	
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	165	125-205	
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	193	147-239	
(CA 15-3)糖类抗原 15-3				
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	63	47.9-78.1	
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	61.2	46.5-75.9	
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	54	41-67	
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	67.5	51.3-83.7	
(CA 19-9)糖类抗原 19-9				
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	366	278-454	
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	166	126-206	
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	82	62.3-102	
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	176	134-218	
(CA 72-4)糖类抗原 72-4				
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	14.4	10.9-17.9	
(CEA)癌胚抗原				
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	44	33.4-54.6	
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	32.6	24.8-40.4	
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	26.4	20.1-32.7	
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	40.5	30.8-50.2	
(Cyfra 21-1)细胞角蛋白19片段				
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	14.3	10.9-17.7	
(Ferritin)铁蛋白				
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	µg/L	404	307-501	
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	µg/L	207	157-257	
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/L	310	236-384	
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence Immunoassay (CLIA)	µg/L	266	202-330	
(F-PSA)游离前列腺特异性抗原				
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.8	3.65-5.95	
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.65	4.29-7.01	

REF 510-2	6x3mL	EXP 2026-07-22	LOT 30A891C	更新日期： 2021-09 2 / 2
------------------	-------	-----------------------	--------------------	------------------------

批号：30A891C				
项目\仪器\试剂方法	单位	均值	+ / - 2 SD	
Roche Modular E170 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.88	2.95-4.81	
Siemens ADVIA Centaur XP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.39	3.34-5.44	
(PSA)前列腺特异性抗原				
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16.3	12.4-20.2	
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	18.2	13.8-22.6	
Roche Modular E170 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	16.6	12.6-20.6	
Siemens ADVIA Centaur XP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16.4	12.5-20.3	
(SCC)鳞状细胞癌抗原				
Abbott AxSYM Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	14.2	10.8-17.6	
(TPA)组织多肽抗原				
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/L	7340	5570-9100	
(TPS)组织多肽特异抗原				
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/L	75	57-93	
(β-hCG)人绒毛膜促性腺激素β亚基				
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	154	117-191	
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	90	68.4-112	

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。