

REF 500-3	6x3mL	EXP 2025-05-01	LOT 2FB6E4C	更新日期: 2024-06
				1 / 3

本报告含以下项目的示范实验室数据

分析项目(Analyte)	分析项目(Analyte)	分析项目(Analyte)
(25-OH-Vit D)25羟基维生素D	(Cortisol)皮质醇	(C-P)C-肽
(E2)雌二醇	(FA)叶酸	(FSH)促卵泡刺激素
(FT3)游离T3	(FT4)游离T4	(hCG)人绒毛膜促性腺激素
(INS)胰岛素	(LH)促黄体生成素	(PRL)催/ 泌乳素
(PROG)孕酮	(TES)睾酮	(TSH)促甲状腺素
(TT3)总三碘甲状原氨酸	(TT4)总T4	(VB12)维生素B12

注: 此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同, 或因制造商检测方法的改变, 均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求, 实验室须遵循相关技术规范确立自己的均值和可接受范围。

REF 500-3	6x3mL	EXP 2025-05-01	LOT 2FB6E4C	更新日期: 2024-06 2 / 3
------------------	-------	-----------------------	--------------------	------------------------

批号: 2FB6E4C				
项目\仪器\试剂方法	单位	均值	+ / - 2 SD	
(25-OH-Vit D)25羟基维生素D				
Abbott ARCHITECT i series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	34.6	26.3-42.9	
Abbott ARCHITECT i series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	nmol/L	86.6	65.8-107	
(Cortisol)皮质醇				
Abbott ARCHITECT i series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	251	191-311	
(C-P)C-肽				
Abbott ARCHITECT i series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.03	6.1-9.96	
(E2)雌二醇				
Abbott ARCHITECT i series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1280	971-1580	
Autolumo A2000/A2000 PLUS				
Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	3230	2450-4000	
Beckman Access				
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	229	183-275	
(FA)叶酸				
Abbott ARCHITECT i series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.97	2.26-3.68	
(FSH)促卵泡刺激素				
Abbott ARCHITECT i series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	37.4	28.4-46.4	
Beckman Access				
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	7.53	6.03-9.03	
Autolumo A2000/A2000 PLUS				
Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/L	42	31.9-52.1	
(FT3)游离T3				
Autolumo A2000/A2000 PLUS				
Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	26.1	19.8-32.4	
Beckman Access				
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	2.4	1.92-2.88	
Abbott ARCHITECT i series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	19.8	15.1-24.5	
(FT4)游离T4				
Abbott ARCHITECT i series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	32.4	24.6-40.2	
Autolumo A2000/A2000 PLUS				
Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	40.5	30.8-50.2	
(hCG)人绒毛膜促性腺激素				
Beckman Access				
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	15.7	12.6-18.8	
Abbott ARCHITECT i series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	508	386-630	
(INS)胰岛素				
Abbott ARCHITECT i series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	99.2	75.4-123	
(LH)促黄体生成素				
Abbott ARCHITECT i series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	18.6	14.1-23.1	
Autolumo A2000/A2000 PLUS				
Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/L	48	36.5-59.5	
Beckman Access				
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	2.05	1.65-2.45	
(PRL)催/泌乳素				
Beckman Access				
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	6.5	5.2-7.8	
Autolumo A2000/A2000 PLUS				
Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	67.1	51-83.2	
Abbott ARCHITECT i series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	45.3	34.4-56.2	
(PROG)孕酮				
Abbott ARCHITECT i series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	22.9	17.4-28.4	

注: 此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同, 或因制造商检测方法的改变, 均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求, 实验室须遵循相关技术规范确立自己的均值和可接受范围。

REF 500-3	6x3mL	EXP 2025-05-01	LOT 2FB6E4C	更新日期： 2024-06 3 / 3
------------------	-------	-----------------------	--------------------	------------------------

Beckman Access Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.91	0.73-1.09
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	26.7	20.3-33.1
(TES)睾酮			
Beckman Access Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.51	0.41-0.61
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12	9.13-14.9
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.79	5.92-9.66
(TSH)促甲状腺素			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	16	12.2-19.8
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	mIU/L	22.9	17.4-28.4
(TT3)总三碘甲状腺原氨酸			
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.01	3.05-4.97
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.99	2.27-3.71
(TT4)总T4			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	nmol/L	135	103-167
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	nmol/L	162	123-201
(VB12)维生素B12			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pg/mL	729	554-904

注：此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的变化，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。