

REF	510-1	6x3mL	EXP	2025-09-26	LOT	30A61AA	更新日期:	2024-08
								1 / 2

		批号: 30A61AA					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数			
(AFP)甲胎蛋白							
Roche cobas E601/602	ng/mL	3.68	9.6%	<5	<100		
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)							
Siemens ADVIA Centaur XP/XPT	ng/mL	2.76	9.8%	<5	<100		
Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)							
Beckman UniCel Dxl 600/800	ng/mL	4.23	9.7%	<5	<100		
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)							
(BNP)B型利钠肽							
AB Sciex API Series	pg/mL	1.82	23.5%	<5	<100		
ABX \ Chemiluminescence enzyme immunoassay (CLEIA)							
(CA 125)糖类抗原 125							
Abbott ARCHITECT i series	IU/mL	32.6	10.1%	<5	<100		
Abbott \ Chemiluminescence Immunoassay (CLIA)							
Siemens ADVIA Centaur XP/XPT	IU/mL	28.3	9.8%	<5	<100		
Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)							
DiaSorin Liaison /XL	IU/mL	32.6	9.7%	<5	<100		
DiaSorin \ Chemiluminescence Immunoassay (CLIA)							
Roche cobas E601/602	IU/mL	29.4	10%	<5	<100		
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)							
(CA 15-3)糖类抗原 15-3							
Abbott ARCHITECT i series	IU/mL	15.3	9.4%	<5	<100		
Abbott \ Chemiluminescence Immunoassay (CLIA)							
DiaSorin Liaison /XL	IU/mL	14.8	9.8%	<5	<100		
DiaSorin \ Chemiluminescence Immunoassay (CLIA)							
Roche cobas E601/602	IU/mL	14.9	9.7%	<5	<100		
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)							
Siemens ADVIA Centaur XP/XPT	IU/mL	14.9	10%	<5	<100		
Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)							
(CA 19-9)糖类抗原 19-9							
Roche cobas E601/602	IU/mL	14.2	9.4%	<5	<100		
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)							
Siemens ADVIA Centaur XP/XPT	IU/mL	26.2	10%	<5	<100		
Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)							
DiaSorin Liaison /XL	IU/mL	27.6	9.1%	<5	<100		
DiaSorin \ Chemiluminescence Immunoassay (CLIA)							
Abbott ARCHITECT i series	IU/mL	45.3	9.7%	<5	<100		
Abbott \ Chemiluminescence Immunoassay (CLIA)							
(CA242)糖类抗原 242							
Manual	U/mL	13.7	9.4%	<5	<100		
GLD \ Enzyme linked immunosorbent assays (ELISA)							
(CEA)癌胚抗原							
Siemens ADVIA Centaur XP/XPT	ng/mL	1.54	9.7%	<5	<100		
Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)							
Roche cobas E601/602	ng/mL	2.73	10.3%	<5	<100		
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)							
Abbott ARCHITECT i series	ng/mL	3.04	10.1%	<5	<100		
Abbott \ Chemiluminescence Immunoassay (CLIA)							
Beckman UniCel Dxl 600/800	ng/mL	2.25	9.3%	<5	<100		
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)							
(Ferritin)铁蛋白							
Siemens ADVIA Centaur XP/XPT	ng/mL	53.7	9.7%	<5	<100		
Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)							
Roche cobas c 501/502/701/702	ng/mL	72.6	9.5%	<5	<100		
Roche \ Immunoturbidimetry (ITA)							
Beckman Access	ng/mL	49.7	9.6%	<5	<100		
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)							
Abbott ARCHITECT i series	ng/mL	62.8	9.9%	<5	<100		
Abbott \ Chemiluminescence Immunoassay (CLIA)							
(F-PSA)游离前列腺特异性抗原							
Beckman Access	ng/mL	1.58	9.9%	<5	<100		
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)							
Siemens ADVIA Centaur XP/XPT	ng/mL	0.94	11%	<5	<100		
Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)							
Roche cobas E601/602	ng/mL	1.11	8.9%	<5	<100		
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)							
Abbott ARCHITECT i series	ng/mL	0.95	8.7%	<5	<100		
Abbott \ Chemiluminescence Immunoassay (CLIA)							
(hCG)人绒毛膜促性腺激素							
Siemens ADVIA Centaur XP/XPT	IU/L	10.6	10.2%	<5	<100		
Siemens \ Enzyme linked immunosorbent assays (ELISA)							
(PSA)前列腺特异性抗原							
Abbott ARCHITECT i series	ng/mL	2.37	9.7%	<5	<100		
Abbott \ Chemiluminescence Immunoassay (CLIA)							

REF	510-1	6x3mL	EXP	2025-09-26	LOT	30A61AA	更新日期:	2024-08
								2 / 2

		批号: 30A61AA			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.28	12%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.24	10.6%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.26	9.9%	<5	<100
(SCC)鳞状细胞癌抗原					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.48	8.2%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.69	9.5%	<5	<100
(TPA)组织多肽抗原					
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	U/L	392	9.9%	<5	<100
(TPS)组织多肽特异抗原					
Manual OTHER \ Enzyme linked immunosorbent assays (ELISA)	U/L	16.4	9.9%	<5	<100
(β-hCG)人绒毛膜促性腺激素β亚基					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.014	9.2%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。