

REF	513-1	6x3mL	EXP	2027-11-26	LOT	30EF654	更新日期:	2024-08
								1 / 4

		批号: 30EF654			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
(AFP)甲胎蛋白					
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.66	3.6%	12	314
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.42	8.8%	<5	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	5.31	8.1%	<5	<100
Beckman UniCel DxI 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.6	9.1%	<5	<100
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.64	10%	<5	<100
Roche Modular E170 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.88	10.6%	<5	<100
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.44	9.9%	<5	<100
(CA 125)糖类抗原 125					
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	39	2.8%	12	266
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	IU/mL	26.3	9.9%	<5	<100
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	22	10%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	30.7	10.5%	<5	<100
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	32.2	8.7%	<5	<100
Roche Modular E170 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	26.5	10.6%	<5	<100
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	24.9	9.6%	<5	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	26.7	8.3%	<5	<100
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	34.5	8.1%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	39.8	10.3%	<5	<100
(CA 15-3)糖类抗原 15-3					
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	14.5	3.6%	12	323
Roche Modular E170 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	13.4	12.4%	<5	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	17.4	7.5%	<5	<100
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	IU/mL	23.1	9.9%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	15.6	11.8%	<5	<100
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	18.5	7%	<5	<100
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	13.4	10%	<5	<100
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	16.6	7.2%	<5	<100
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	15.4	7.2%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	15.2	11.1%	<5	<100
(CA 19-9)糖类抗原 19-9					
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	43.8	6.4%	12	292
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	27.4	8%	<5	<100
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	15.5	7.1%	<5	<100
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	47.9	6.9%	<5	<100
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	23.3	9.4%	<5	<100
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	IU/mL	20.4	9.9%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	23.5	11.7%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	44.3	9.2%	<5	<100
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	14.6	10%	<5	<100
Roche Modular E170 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	11.6	11.2%	<5	<100

REF	513-1	6x3mL	EXP	2027-11-26	LOT	30EF654	更新日期:	2024-08
								2 / 4

		批号: 30EF654					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数			
(CA 50)糖类抗原 50							
Mindray CL-8000i							
Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	10.4	10%	<5	<100		
HYBIOME EVERESYS Series							
HYBIOME \ Chemiluminescence Immunoassay (CLIA)	IU/mL	17.3	9.9%	<5	<100		
(CA 72-4)糖类抗原 72-4							
Roche Modular E170							
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	2.1	8.1%	<5	<100		
Roche Modular E170							
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	3.57	12.3%	<5	<100		
Maccura i3000							
MAKER \ Chemiluminescence Immunoassay (CLIA)	IU/mL	3.58	4.7%	<5	<100		
HYBIOME EVERESYS Series							
HYBIOME \ Chemiluminescence Immunoassay (CLIA)	IU/mL	1.05	9.9%	<5	<100		
(CA242)糖类抗原 242							
HYBIOME EVERESYS Series							
HYBIOME \ Chemiluminescence Immunoassay (CLIA)	U/mL	9.51	9.9%	<5	<100		
Maccura i3000							
MAKER \ Chemiluminescence Immunoassay (CLIA)	U/mL	9.6	1.5%	<5	<100		
Mindray CL-8000i							
Mindray \ Chemiluminescence Immunoassay (CLIA)	U/mL	7.64	10%	<5	<100		
Siemens Immulite1000/2000							
Siemens \ Chemiluminescence Immunoassay (CLIA)	U/mL	9.89	9.9%	<5	<100		
(CEA)癌胚抗原							
Abbott Alinity i Series							
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.94	7.3%	10	257		
Abbott ARCHITECT i series							
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.23	11.6%	<5	<100		
Siemens ADVIA Centaur XP/XPT							
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.5	7.2%	<5	<100		
HYBIOME EVERESYS Series							
HYBIOME \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.61	9.9%	<5	<100		
Mindray CL-8000i							
Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.96	10%	<5	<100		
Roche Modular E170							
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.75	12.2%	<5	<100		
Siemens ADVIA Centaur XP/XPT							
Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.54	12.3%	<5	<100		
Abbott ARCHITECT i series							
Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.99	8.7%	<5	<100		
Roche Modular E170							
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.42	7.6%	<5	<100		
Beckman UniCel DxI 600/800							
Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.54	8.3%	<5	<100		
(Cyfra 21-1)细胞角蛋白19片段							
Abbott Alinity i Series							
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.47	6.1%	12	327		
Maccura i3000							
MAKER \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.1	1.7%	<5	<100		
Roche Modular E170							
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.94	5.7%	<5	<100		
Mindray CL-8000i							
Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.33	10%	<5	<100		
Roche Modular E170							
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.48	12.1%	<5	<100		
(Ferritin)铁蛋白							
Abbott ARCHITECT i series							
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	88.7	12.3%	<5	<100		
Roche Modular E170							
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	53.4	9.2%	<5	<100		
Abbott ARCHITECT i series							
Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	48.3	8.9%	<5	<100		
Roche Modular E170							
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	94.2	12.2%	<5	<100		
Beckman Access							
Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	35.7	9%	<5	<100		
HYBIOME EVERESYS Series							
HYBIOME \ Chemiluminescence Immunoassay (CLIA)	ng/mL	43.8	9.9%	<5	<100		
Mindray CL-8000i							
Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	46	10%	<5	<100		
Siemens ADVIA Centaur XP/XPT							
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	45.6	9.4%	<5	<100		
(F-PSA)游离前列腺特异性抗原							
Abbott Alinity i Series							
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.691	3.3%	12	328		
Siemens ADVIA Centaur XP/XPT							
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.15	7.8%	<5	<100		
HYBIOME EVERESYS Series							
HYBIOME \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.929	9.9%	<5	<100		
Abbott ARCHITECT i series							
Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.949	9.5%	<5	<100		

REF	513-1	6x3mL	EXP	2027-11-26	LOT	30EF654	更新日期:	2024-08
								3 / 4

		批号: 30EF654			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
Roche Modular E170 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.67	10.9%	<5	<100
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.93	10%	<5	<100
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.02	8.8%	<5	<100
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.42	6.3%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.82	10.2%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.67	12.5%	<5	<100
(hCG)人绒毛膜促性腺激素					
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	19.4	8.3%	<5	<100
HYBIO ME EVERESYS Series HYBIO ME \ Chemiluminescence Immunoassay (CLIA)	IU/L	6.09	9.9%	<5	<100
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/L	8.73	10%	<5	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Enzyme linked immunosorbent assays (ELISA)	IU/L	14.4	9.7%	<5	<100
Roche Modular E170 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	15.6	11.7%	<5	<100
(HE4)人附睾蛋白4					
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	44.6	6.4%	11	283
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	36.8	7.6%	<5	<100
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	38.2	10%	<5	<100
Roche Modular E170 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	45.9	12.5%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	pmol/L	46.2	5.1%	<5	<100
(PSA)前列腺特异性抗原					
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.21	5%	12	323
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.48	8.9%	<5	<100
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.43	9%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.78	11.8%	<5	<100
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.5	8.8%	<5	<100
Roche Modular E170 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.69	12.3%	<5	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.3	9.6%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.48	12.2%	<5	<100
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.83	10%	<5	<100
HYBIO ME EVERESYS Series HYBIO ME \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.84	9.9%	<5	<100
(SCC)鳞状细胞癌抗原					
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.4	5.5%	12	317
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.24	8.9%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.37	3.1%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.53	11.8%	<5	<100
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.25	10%	<5	<100
(TPA)组织多肽抗原					
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/L	494	9.7%	<5	<100
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	U/L	278	10.4%	<5	<100
(TPS)组织多肽特异抗原					
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	U/L	16.4	12.4%	<5	<100
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/L	36.7	8.7%	<5	<100
(β-hCG)人绒毛膜促性腺激素β亚基					
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.033	9.3%	<5	<100

REF	513-1	6x3mL	EXP	2027-11-26	LOT	30EF654	更新日期：	2024-08
								4 / 4

		批号： 30EF654	
项目\仪器\试剂\方法	单位	均值 \ CV	Lab \ 点数

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。