

REF 521-1	6x3mL	EXP 2024-04-04	LOT 31B2AE8	更新日期: 2024-08
				1 / 4

批号: 31B2AE8					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
(Alb)白蛋白					
Beckman AU series Beckman Coulter \ Bromocresol purple	g/L	22.4	11.1%	<5	<100
(ASO)抗链球菌溶血素 O 抗体					
Beckman AU series BSBE \ Immunoturbidimetry (ITA)	IU/mL	127	3.8%	9	256
Beckman Immage 800 Beckman Coulter \ OTHER	IU/mL	90.2	5.4%	8	175
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	151	11.9%	6	<100
TOSHIBA TBA OTHER \ OTHER	IU/mL	150	7%	<5	<100
Siemens BN ProSpec Dedicated \ Immunonephelometry (INA)	IU/mL	152	2.9%	21	313
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	IU/mL	251	4.9%	<5	<100
Mindray BS series Whitman \ Immunoturbidimetry (ITA)	IU/mL	270	1%	<5	<100
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	IU/mL	158	5.2%	16	648
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	IU/mL	96.6	3.8%	13	<100
TOSHIBA TBA MedicalSystem \ Immunoturbidimetry (ITA)	IU/mL	165	2.9%	12	332
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	155	3.2%	<5	<100
(AT-III)抗凝血酶III					
Beckman Immage 800 Beckman Coulter \ Nephelometry	g/L	0.136	11.3%	<5	<100
(C1q)补体C1q					
Beckman AU series chongqiong zybio \ OTHER	mg/L	133	8.5%	11	182
(C3)补体 C3					
Beckman AU series BSBE \ Immunoturbidimetry (ITA)	mg/dL	66.9	6.9%	8	243
Beckman Immage 800 Beckman Coulter \ OTHER	mg/dL	64.4	5.5%	8	191
Siemens BN II Dedicated \ Immunonephelometry (INA)	mg/dL	71.1	3.6%	5	118
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	mg/dL	70	5.2%	<5	<100
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	mg/dL	65.2	3.1%	17	<100
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	mg/dL	70.5	2.8%	16	664
Siemens BN ProSpec Dedicated \ Immunonephelometry (INA)	mg/dL	75.6	1.7%	16	196
TOSHIBA TBA LEADMAN \ Immunoturbidimetry (ITA)	mg/dL	78.7	5.3%	13	382
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	76.7	11.7%	<5	<100
(C4)补体 C4					
Beckman AU series BSBE \ Immunoturbidimetry (ITA)	mg/dL	16.8	3.7%	9	252
Beckman Immage 800 Beckman Coulter \ OTHER	mg/dL	18.5	5.5%	8	187
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	mg/dL	12.8	5.5%	5	<100
Siemens BN II Dedicated \ Immunonephelometry (INA)	mg/dL	17.1	2.7%	5	118
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	mg/dL	27.4	5.8%	<5	<100
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	mg/dL	17.8	4%	17	<100
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	mg/dL	15.2	2.8%	16	666
Siemens BN ProSpec Dedicated \ Immunonephelometry (INA)	mg/dL	18.3	3.4%	16	196
TOSHIBA TBA LEADMAN \ Immunoturbidimetry (ITA)	mg/dL	16.1	4.2%	13	357
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	19.5	11.4%	<5	<100
(Cer)铜蓝蛋白					
Siemens BN II Dedicated \ Immunonephelometry (INA)	g/L	0.185	3.3%	5	118
Siemens BN ProSpec Dedicated \ Immunonephelometry (INA)	g/L	0.192	2.2%	13	<100
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	g/L	0.146	6.4%	10	<100
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	g/L	0.263	10.9%	<5	<100
(CRP)C反应蛋白					
Beckman Immage 800 Beckman Coulter \ OTHER	mg/dL	1.28	4.2%	8	190
Beckman AU series DongOu \ Immunoturbidimetry (ITA)	mg/dL	1.37	2.1%	5	125

REF	521-1	6x3mL	EXP	2024-04-04	LOT	31B2AE8	更新日期:	2024-08
								2 / 4

		批号: 31B2AE8			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	1.44	11.2%	<5	<100
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	mg/dL	1.33	4%	17	<100
Siemens BN ProSpec Dedicated \ Immunonephelometry (INA)	mg/dL	1.32	3.1%	16	196
Mindray BC-5100/BC-5300 Series Dedicated \ Immunoturbidimetry (ITA)	mg/dL	1.43	2.9%	15	249
TOSHIBA TBA BSBE \ Immunoturbidimetry (ITA)	mg/dL	1.39	5.7%	12	337
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	1.46	4.9%	<5	<100
(Ferritin)铁蛋白					
Siemens BN II Dedicated \ Immunoturbidimetry (ITA)	ng/mL	48.2	12.7%	<5	<100
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	ng/mL	54.3	12.6%	<5	<100
(Free Kappa Light Chain)游离轻链κ					
Siemens BN II Dedicated \ Immunonephelometry (INA)	mg/L	14	4.4%	5	118
(Free Lambda Light Chain)游离轻链λ					
Siemens BN II Dedicated \ Immunonephelometry (INA)	mg/L	17.2	4.1%	5	118
(Hap (HPT))触珠蛋白 (结合珠蛋白)					
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	mg/L	695	10.5%	<5	<100
(hs-CRP)超敏C反应蛋白					
TOSHIBA TBA OTHER \ OTHER	mg/L	16.8	2.2%	5	<100
(IgA)免疫球蛋白 IgA					
Beckman Immage 800 Beckman Coulter \ OTHER	mg/dL	109	4.7%	8	191
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	mg/dL	119	4%	5	<100
Siemens BN II Dedicated \ Immunonephelometry (INA)	mg/dL	125	3.2%	5	118
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	mg/dL	210	4.2%	<5	<100
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	mg/dL	110	2.9%	17	<100
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	mg/dL	121	3.4%	16	664
Siemens BN ProSpec Dedicated \ Immunonephelometry (INA)	mg/dL	130	2.9%	16	196
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	148	10.8%	<5	<100
(IgE)免疫球蛋白 IgE					
Siemens BN II Dedicated \ Immunonephelometry (INA)	IU/mL	28.2	4.2%	5	118
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	43.5	11%	<5	<100
(IgG)免疫球蛋白 IgG					
Beckman Immage 800 Beckman Coulter \ OTHER	mg/dL	671	5.2%	9	191
Siemens BN II Dedicated \ Immunonephelometry (INA)	mg/dL	758	2.9%	5	118
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	mg/dL	590	3%	<5	<100
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	mg/dL	878	3.3%	<5	<100
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	mg/dL	688	3.8%	17	<100
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	mg/dL	685	2.1%	16	666
Siemens BN ProSpec Dedicated \ Immunonephelometry (INA)	mg/dL	741	3.9%	16	196
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	895	10.7%	<5	<100
(IgG4)免疫球蛋白G4亚型					
Siemens BN II Dedicated \ Immunonephelometry (INA)	g/L	363	2.5%	5	118
(IgM)免疫球蛋白 IgM					
Beckman Immage 800 Beckman Coulter \ OTHER	mg/dL	84.3	5.8%	8	189
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	mg/dL	71.3	7.2%	5	<100
Siemens BN II Dedicated \ Immunonephelometry (INA)	mg/dL	87.1	2.9%	5	118
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	mg/dL	67.2	4.8%	<5	<100
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	mg/dL	86.2	3.4%	17	<100
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	mg/dL	78.5	2.6%	16	658

REF	521-1	6x3mL	EXP	2024-04-04	LOT	31B2AE8	更新日期:	2024-08
								3 / 4

		批号: 31B2AE8					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数			
Siemens BN ProSpec Dedicated \ Immunonephelometry (INA)	mg/dL	86.5	2.9%	16	196		
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	103	11.2%	<5	<100		
(Kappa Light Chain)κ轻链							
Siemens BN II Dedicated \ Immunonephelometry (INA)	g/L	1.75	4.8%	5	118		
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	g/L	5.44	9.5%	17	<100		
Siemens BN ProSpec Dedicated \ Immunonephelometry (INA)	g/L	1.76	2.2%	16	196		
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	g/L	7.06	10.4%	<5	<100		
(Lambda Light Chain)λ 轻链							
Siemens BN II Dedicated \ Immunonephelometry (INA)	g/L	0.885	4.9%	5	119		
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	g/L	2.86	8.3%	17	<100		
Siemens BN ProSpec Dedicated \ Immunonephelometry (INA)	g/L	1.01	2.6%	16	197		
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	g/L	3.66	10.6%	<5	<100		
(PAB)前白蛋白							
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	g/L	0.155	2.2%	<5	<100		
Beckman AU series DongOu \ Immunoturbidimetry (ITA)	g/L	0.155	2%	14	298		
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	g/L	0.169	11.8%	<5	<100		
(PRO/ TP)(总) 蛋白							
Beckman AU series Beckman Coulter \ OTHER	g/L	6.97	2.9%	11	183		
(RBP)视黄醇结合蛋白							
Siemens ADVIA 1800/2400 BIOSINO \ Immunoturbidimetry (ITA)	mg/dL	2.81	4.2%	16	653		
Beckman AU series DongOu \ Immunoturbidimetry (ITA)	mg/dL	2.63	2.2%	14	297		
Beckman AU series DongOu \ OTHER	mg/dL	2.5	6%	11	183		
Siemens BN II Dedicated \ Immunoturbidimetry (ITA)	mg/dL	2.56	12.4%	<5	<100		
(RF)类风湿因子							
Beckman AU series BSBE \ Immunoturbidimetry (ITA)	IU/mL	39.4	3.4%	9	252		
Beckman Immage 800 Beckman Coulter \ OTHER	IU/mL	34.7	5.4%	8	143		
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	37.8	2.5%	5	<100		
TOSHIBA TBA OTHER \ OTHER	IU/mL	48.4	2.9%	<5	<100		
Siemens BN ProSpec Dedicated \ Immunonephelometry (INA)	IU/mL	33.5	6%	21	313		
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	IU/mL	36.7	6.5%	16	662		
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	IU/mL	27.9	7%	13	<100		
TOSHIBA TBA BSBE \ Immunoturbidimetry (ITA)	IU/mL	40.1	2.7%	12	326		
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	40.6	1.2%	<5	<100		
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	32.6	11.6%	<5	<100		
(TRF)转铁蛋白							
Beckman Immage 800 Beckman Coulter \ OTHER	mg/dL	124	6.7%	9	155		
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	mg/dL	121	2.7%	8	327		
Siemens BN ProSpec Dedicated \ Immunonephelometry (INA)	mg/dL	122	3%	16	197		
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	116	12.5%	<5	<100		
(α-1-AG)α-1-酸性糖蛋白							
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	g/L	0.432	12.4%	<5	<100		
(α-1-AT)α-1-抗胰蛋白酶							
Beckman AU series OTHER \ Immunoturbidimetry (ITA)	mg/dL	75.3	4.8%	14	297		
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	72.5	11.3%	<5	<100		
(α-1-G)α-1-球蛋白							
TOSHIBA TBA LEADMAN \ Capillary Method	%	10.1	6.4%	7	212		

REF	521-1	6x3mL	EXP	2024-04-04	LOT	31B2AE8	更新日期：	2024-08 4 / 4
------------	-------	-------	------------	------------	------------	---------	-------	------------------

		批号： 31B2AE8			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(α2M)α-2-巨球蛋白					
Beckman Image 800 Beckman Coulter \ Immunoturbidimetry (ITA)	mg/L	1070	11.2%	<5	<100
(β-2-G)β-2-球蛋白					
TOSHIBA TBA LEADMAN \ Capillary Method	%	2.02	6.4%	<5	<100
(β2M)β-2-微球蛋白					
Beckman AU series Dedicated \ Immunoturbidimetry (ITA)	ng/mL	1613	5.8%	9	255
Roche Hitachi/Hitachi Systems OTHER \ OTHER	ng/mL	1771	18.1%	5	<100
TOSHIBA TBA OTHER \ OTHER	ng/mL	2133	7.5%	<5	<100
Beckman AU series LEADMAN \ Immunoturbidimetry (ITA)	ng/mL	2389	2.7%	14	300
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	ng/mL	2517	7.1%	12	337
Roche Hitachi/Hitachi Systems MedicalSystem \ Immunoturbidimetry (ITA)	ng/mL	4466	7.4%	11	231
Beckman AU series DongOu \ OTHER	ng/mL	2220	6.3%	11	182
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	ng/mL	2800	12.3%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。