

REF	521-1	6x3mL	EXP	2026-05-22	LOT	31B2BB1	更新日期:	2024-08
								1 / 3

		批号: 31B2BB1			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
(Alb)白蛋白					
Olympus AU400					
Olympus \ Bromocresol green	g/L	20.6	10%	<5	<100
(ASO)抗链球菌溶血素 O 抗体					
Mindray BS series					
MedicalSystem \ OTHER	IU/mL	146	6.1%	6	135
Siemens BN ProSpec					
Dedicated \ Immunonephelometry (INA)	IU/mL	131	2.6%	<5	<100
Beckman Immage 800					
Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	86.7	4.5%	<5	<100
Siemens ADVIA 1800/2400					
DiaSys \ Immunoturbidimetry (ITA)	IU/mL	142	3.7%	<5	135
Beckman AU series					
chongqiong zybio \ Immunoturbidimetry (ITA)	IU/mL	117	2.7%	12	162
Olympus AU400					
Olympus \ Immunonephelometry (INA)	IU/mL	151	10%	<5	<100
Siemens ADVIA Centaur XP/XPT					
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	124	0.8%	<5	<100
Beckman AU series					
chongqiong zybio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	120	10%	<5	<100
(AT-III)抗凝血酶III					
Beckman Immage 800					
Beckman Coulter \ Nephelometry	g/L	0.128	10%	<5	<100
(C3)补体 C3					
Beckman AU series					
chongqiong zybio \ Immunoturbidimetry (ITA)	mg/dL	77.3	9.7%	7	130
Siemens ADVIA 1800/2400					
DiaSys \ Immunoturbidimetry (ITA)	mg/dL	76.7	5.2%	6	270
Mindray BS series					
MedicalSystem \ OTHER	mg/dL	67.5	3.4%	6	135
Siemens BN II					
Dedicated \ Immunonephelometry (INA)	mg/dL	81	3.3%	<5	<100
Beckman Immage 800					
Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	72	4.5%	<5	<100
Olympus AU400					
Olympus \ Immunonephelometry (INA)	mg/dL	72.4	10%	<5	<100
Siemens ADVIA Centaur XP/XPT					
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	mg/dL	0.738	2%	<5	<100
(C4)补体 C4					
Beckman AU series					
chongqiong zybio \ Immunoturbidimetry (ITA)	mg/dL	18.4	7.3%	7	112
Siemens ADVIA 1800/2400					
DiaSys \ Immunoturbidimetry (ITA)	mg/dL	17	6%	6	272
Mindray BS series					
MedicalSystem \ OTHER	mg/dL	17.3	1.9%	6	134
Siemens BN II					
Dedicated \ Immunonephelometry (INA)	mg/dL	19.7	3.8%	<5	<100
Beckman Immage 800					
Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	19.5	5.4%	<5	<100
Siemens ADVIA Centaur XP/XPT					
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	mg/dL	0.162	4.6%	<5	<100
Olympus AU400					
Olympus \ Immunonephelometry (INA)	mg/dL	19.3	10%	<5	<100
(Cer)铜蓝蛋白					
Siemens BN ProSpec					
Dedicated \ Immunonephelometry (INA)	g/L	0.173	1.5%	<5	<100
Olympus AU400					
Olympus \ Immunonephelometry (INA)	g/L	0.215	10%	<5	<100
(CRP)C反应蛋白					
Mindray BS series					
MedicalSystem \ OTHER	mg/dL	1.44	15.1%	6	135
Olympus AU400					
Olympus \ Immunonephelometry (INA)	mg/dL	1.53	10%	<5	<100
(Ferritin)铁蛋白					
Siemens BN II					
Siemens \ Immunonephelometry (INA)	ng/mL	49.8	10%	<5	<100
Olympus AU400					
Olympus \ Immunoturbidimetry (ITA)	ng/mL	51	10%	<5	<100
(Hap (HPT))触珠蛋白 (结合珠蛋白)					
Olympus AU400					
Olympus \ Enzyme linked immunosorbent assays (ELISA)	mg/L	760	10%	<5	<100
(IgA)免疫球蛋白 IgA					
Siemens ADVIA 1800/2400					
DiaSys \ Immunoturbidimetry (ITA)	mg/dL	129	8%	6	270
Siemens BN II					
Dedicated \ Immunonephelometry (INA)	mg/dL	141	3.3%	<5	<100
Beckman Immage 800					
Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	115	4.7%	<5	<100
Olympus AU400					
Olympus \ Immunonephelometry (INA)	mg/dL	139	10%	<5	<100
Siemens ADVIA Centaur XP/XPT					
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	mg/dL	121	4.2%	<5	<100

REF	521-1	6x3mL	EXP	2026-05-22	LOT	31B2BB1	更新日期:	2024-08
								2 / 3

批号: 31B2BB1						
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数		
(IgE)免疫球蛋白 IgE						
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	IU/mL	21.9	10%	<5	<100	
(IgG)免疫球蛋白 IgG						
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	mg/dL	677	7%	6	273	
Siemens BN II Dedicated \ Immunonephelometry (INA)	mg/dL	750	4%	<5	<100	
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	684	4.4%	<5	<100	
Siemens ADVIA Centaur XP/XPT Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	mg/dL	775	3.6%	<5	<100	
Olympus AU400 Olympus \ Immunonephelometry (INA)	mg/dL	729	10%	<5	<100	
(IgM)免疫球蛋白 IgM						
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	mg/dL	82.3	6.4%	6	269	
Siemens BN II Dedicated \ Immunonephelometry (INA)	mg/dL	90.5	2.9%	<5	<100	
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	82.6	4.9%	<5	<100	
Siemens ADVIA Centaur XP/XPT Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	mg/dL	72.4	3.3%	<5	<100	
Olympus AU400 Olympus \ Immunonephelometry (INA)	mg/dL	76	10%	<5	<100	
(Kappa Light Chain)κ轻链						
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	g/L	5.75	10%	<5	<100	
(Lambda Light Chain)λ 轻链						
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	g/L	3.24	10%	<5	<100	
(PAB)前白蛋白						
Siemens ADVIA 1800/2400 BIOSINO \ Immunoturbidimetry (ITA)	g/L	182	5.9%	6	268	
Mindray BS series OTHER \ OTHER	g/L	0.175	2.4%	6	135	
Beckman AU series chongqiong zybio \ Immunoturbidimetry (ITA)	g/L	0.018	9.5%	11	145	
Siemens ADVIA Centaur XP/XPT Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	g/L	172	0.8%	<5	<100	
Olympus AU400 Olympus \ Immunonephelometry (INA)	g/L	0.204	10%	<5	<100	
(RBP)视黄醇结合蛋白						
Siemens ADVIA 1800/2400 BIOSINO \ Immunoturbidimetry (ITA)	mg/dL	3.06	5.1%	6	274	
Siemens BN II Siemens \ Immunonephelometry (INA)	mg/dL	3.16	10%	<5	<100	
(RF)类风湿因子						
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	IU/mL	37.9	7%	6	270	
Mindray BS series MedicalSystem \ OTHER	IU/mL	42.5	2%	6	118	
Siemens BN ProSpec Dedicated \ Immunonephelometry (INA)	IU/mL	38.6	3.5%	<5	<100	
Beckman AU series chongqiong zybio \ Immunoturbidimetry (ITA)	IU/mL	39.8	7.4%	13	203	
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	IU/mL	46	10%	<5	<100	
(TRF)转铁蛋白						
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	mg/dL	120	1.4%	<5	136	
Olympus AU400 Olympus \ Immunonephelometry (INA)	mg/dL	133	10%	<5	<100	
(α-1-AG)α-1-酸性糖蛋白						
Olympus AU400 Olympus \ Enzyme linked immunosorbent assays (ELISA)	g/L	0.455	10%	<5	<100	
(α-1-AT)α-1-抗胰蛋白酶						
Olympus AU400 Olympus \ Immunoturbidimetry (ITA)	mg/dL	79.2	10%	<5	<100	
(α2M)α-2-巨球蛋白						
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	mg/L	980	10%	<5	<100	
(β2M)β-2-微球蛋白						
Mindray BS series MedicalSystem \ OTHER	ng/mL	2503	1.6%	6	135	
Beckman AU series chongqiong zybio \ Immunoturbidimetry (ITA)	ng/mL	2080	7.7%	13	212	
Olympus AU400 Olympus \ Immunonephelometry (INA)	ng/mL	2500	12%	<5	<100	

报告类别： 汇总实验室数据报告



质控品名： 特定蛋白液体非定值质控品

REF	521-1	6x3mL	EXP	2026-05-22	LOT	31B2BB1	更新日期：	2024-08
								3 / 3

		批号：31B2BB1	
项目\仪器\试剂\方法	单位	均值 \ CV	Lab \ 点数

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。