

<b>REF</b>	521-2	6x3mL	<b>EXP</b>	2024-04-04	<b>LOT</b>	31B51F8	更新日期:	2024-08
								1 / 6

批号: 31B51F8						
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数		
<b>(Alb)白蛋白</b>						
Siemens ADVIA Centaur CP Dedicated \ Bromocresol purple	g/L	33.1	1.7%	<5	<100	
Roche Hitachi/Hitachi Systems Roche \ Bromocresol purple	g/L	33.2	11.1%	<5	<100	
<b>(ASO)抗链球菌溶血素 O 抗体</b>						
Roche Hitachi/Hitachi Systems BSBE \ Immunoturbidimetry (ITA)	IU/mL	224	2%	9	227	
TOSHIBA TBA MedicalSystem \ Immunoturbidimetry (ITA)	IU/mL	271	2.3%	9	173	
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	IU/mL	217	2.4%	5	145	
Roche cobas c 501/502/701/702 Dedicated \ Immunoturbidimetry (ITA)	IU/mL	276	3.4%	5	144	
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	IU/mL	148	8.9%	46	987	
Beckman Image 800 OTHER \ Immunonephelometry (INA)	IU/mL	131	5.3%	41	1225	
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	IU/mL	234	7.2%	28	1108	
Mindray BS series Mindray \ Immunonephelometry (INA)	IU/mL	261	11.5%	22	374	
Beckman Image 800 Dedicated \ Immunonephelometry (INA)	IU/mL	136	7.9%	21	550	
Beckman Image 800 Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	146	20.7%	20	436	
Beckman AU series HuaChenBio \ Immunonephelometry (INA)	IU/mL	216	1.5%	<5	<100	
Siemens ADVIA Centaur CP Dedicated \ Immunoturbidimetry (ITA)	IU/mL	253	1.6%	<5	<100	
Mindray BS series Whitman \ Immunoturbidimetry (ITA)	IU/mL	165	0.9%	<5	<100	
Siemens ADVIA 1800/2400 Siemens \ Immunoturbidimetry (ITA)	IU/mL	252	1.9%	<5	133	
TOSHIBA TBA MedicalSystem \ Immunonephelometry (INA)	IU/mL	237	8.5%	19	251	
Siemens BN II Siemens \ OTHER	IU/mL	241	6.4%	16	419	
Beckman Image 800 Beckman Coulter \ OTHER	IU/mL	126	12.4%	14	307	
Abbott ARCHITECT c series BSBE \ OTHER	IU/mL	208	11.7%	11	213	
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	239	11.4%	10	225	
<b>(C3)补体 C3</b>						
Roche Hitachi/Hitachi Systems BSBE \ Immunoturbidimetry (ITA)	mg/dL	150	3.8%	9	221	
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	144	3.8%	9	208	
Roche cobas c 501/502/701/702 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	151	3.3%	5	142	
Beckman Image 800 OTHER \ Immunonephelometry (INA)	mg/dL	131	4.2%	43	1257	
Beckman Image 800 Dedicated \ Immunonephelometry (INA)	mg/dL	131	5.3%	<5	<100	
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	135	4%	34	676	
Abbott ARCHITECT c series BSBE \ OTHER	mg/dL	141	1.3%	<5	<100	
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	mg/dL	138	3.8%	28	1121	
Biostec A15/A25 Biostec \ Immunonephelometry (INA)	mg/dL	137	2.6%	27	716	
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	156	6.2%	26	384	
Beckman Image 800 Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	130	5.3%	20	462	
Siemens ADVIA 1800/2400 Siemens \ Immunoturbidimetry (ITA)	mg/dL	135	2.3%	<5	<100	
Siemens ADVIA Centaur CP Dedicated \ Immunoturbidimetry (ITA)	mg/dL	140	2.4%	<5	<100	
Roche Hitachi/Hitachi Systems YongHe SUN \ Immunoturbidimetry (ITA)	mg/dL	130	3.1%	<5	<100	
TOSHIBA TBA MedicalSystem \ Immunoturbidimetry (ITA)	mg/dL	129	10.1%	19	244	
Beckman Image 800 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	139	4.2%	17	495	
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	mg/dL	132	5.1%	16	471	
Siemens BN II Siemens \ OTHER	mg/dL	158	8.9%	16	438	
TOSHIBA TBA LEADMAN \ Immunoturbidimetry (ITA)	mg/dL	156	6.9%	16	342	
<b>(C4)补体 C4</b>						
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	39.2	8.4%	9	214	
Roche Hitachi/Hitachi Systems BSBE \ Immunoturbidimetry (ITA)	mg/dL	38.5	2.1%	9	210	
Roche cobas c 501/502/701/702 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	31.2	4.1%	5	143	

<b>REF</b>	521-2	6x3mL	<b>EXP</b>	2024-04-04	<b>LOT</b>	31B51F8	更新日期:	2024-08
								2 / 6

		批号: 31B51F8					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数			
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	mg/dL	37.2	6%	43	1255		
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	mg/dL	37.5	6.2%	<5	<100		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	35.9	4.9%	34	678		
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	mg/dL	32.3	3.3%	28	1125		
Biostec A15/A25 Biostec \ Immunonephelometry (INA)	mg/dL	42.2	5.1%	27	717		
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	38.5	3.2%	26	404		
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	35.9	5.9%	20	456		
Siemens ADVIA 1800/2400 Siemens \ Immunoturbidimetry (ITA)	mg/dL	36.2	7.3%	<5	<100		
Siemens ADVIA Centaur CP Dedicated \ Immunoturbidimetry (ITA)	mg/dL	40.3	1.9%	<5	<100		
Abbott ARCHITECT c series BSBE \ OTHER	mg/dL	30.8	1.9%	<5	<100		
TOSHIBA TBA MedicalSystem \ Immunoturbidimetry (ITA)	mg/dL	34.2	9.7%	18	255		
Beckman Immage 800 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	37.6	3.9%	17	496		
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	mg/dL	40.9	7.7%	16	470		
Siemens BN II Siemens \ OTHER	mg/dL	35.4	7.7%	16	435		
Beckman Immage 800 Beckman Coulter \ OTHER	mg/dL	37.7	5.7%	14	303		
TOSHIBA TBA LEADMAN \ Immunoturbidimetry (ITA)	mg/dL	34.1	11%	14	299		
<b>(Cer)铜蓝蛋白</b>							
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	g/L	0.337	11%	9	219		
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	g/L	0.302	8%	21	612		
Siemens BN II Siemens \ OTHER	g/L	0.477	8.5%	16	427		
Siemens BN II Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	g/L	0.321	10%	<5	<100		
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	g/L	0.295	6.5%	<5	<100		
Roche cobas 4000/6000/8000 Roche \ Immunonephelometry (INA)	g/L	0.367	0%	<5	<100		
<b>(CRP)C反应蛋白</b>							
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	mg/dL	2.24	3.9%	43	1254		
Roche cobas c 501/502/701/702 Roche \ Immunoturbidimetry (ITA)	mg/dL	2.81	4.3%	<5	<100		
Roche Hitachi/Hitachi Systems chongqiong zybio \ Immunoturbidimetry (ITA)	mg/dL	2.76	7.2%	<5	<100		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	2.45	4.8%	34	676		
Siemens ADVIA 1800/2400 Siemens \ Immunoturbidimetry (ITA)	mg/dL	2.61	3.6%	<5	109		
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	mg/dL	2.39	6.6%	19	525		
Siemens BN II Siemens \ OTHER	mg/dL	2.43	8.7%	16	423		
TOSHIBA TBA BSBE \ Immunoturbidimetry (ITA)	mg/dL	3.07	9.2%	16	333		
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	2.31	5.4%	12	271		
NIHON KOHDEN MEK Series MAKER \ Immunoturbidimetry (ITA)	mg/dL	1.87	2.8%	<5	<100		
Beckman AU series KHB \ Immunoturbidimetry (ITA)	mg/dL	2.08	9.6%	<5	<100		
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	2.84	1.1%	<5	<100		
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mg/dL	2.09	10.5%	<5	<100		
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	2.63	12.2%	<5	<100		
<b>(Ferritin)铁蛋白</b>							
Roche Hitachi/Hitachi Systems MedicalSystem \ Immunoturbidimetry (ITA)	ng/mL	142	3.1%	<5	<100		
Roche Hitachi/Hitachi Systems LEADMAN \ Immunoturbidimetry (ITA)	ng/mL	117	2.7%	17	338		
Siemens BN II Siemens \ OTHER	ng/mL	131	8.8%	16	430		
<b>(Free Kappa Light Chain)游离轻链k</b>							
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	mg/L	10900	9.8%	14	268		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/L	9798	9.6%	12	238		
<b>(Free Lambda Light Chain)游离轻链λ</b>							
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/L	5497	10.5%	26	499		

<b>REF</b>	521-2	6x3mL	<b>EXP</b>	2024-04-04	<b>LOT</b>	31B51F8	更新日期:	2024-08
								3 / 6

		批号: 31B51F8					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数			
<b>(Hap (HPT) )触珠蛋白 (结合珠蛋白)</b>							
Siemens BN II Siemens \ OTHER	mg/L	1432	10.3%	16	437		
Siemens BN II Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mg/L	1369	9.9%	<5	<100		
<b>(HbA1c (%))糖化血红蛋白A1c (%)</b>							
Roche cobas c 501/502/701/702 Dedicated \ immunoturbidimetric	%	5.77	1.3%	5	146		
<b>(hs-CRP)超敏C反应蛋白</b>							
Abbott ARCHITECT c series BSBE \ OTHER	mg/L	29	8.4%	22	455		
Siemens ADVIA Centaur CP OTHER \ Immunoturbidimetry (ITA)	mg/L	27.3	13.5%	<5	<100		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/L	25.7	9.6%	12	274		
Siemens ADVIA Centaur CP Dedicated \ Immunoturbidimetry (ITA)	mg/L	25	2.2%	<5	<100		
Mindray BC-5000 Series Mindray \ Immunonephelometry (INA)	mg/L	27.5	4.1%	<5	<100		
Roche cobas 4000/6000/8000 Roche \ Immunonephelometry (INA)	mg/L	32.1	0%	<5	<100		
<b>(IgA)免疫球蛋白 IgA</b>							
Roche Hitachi/Hitachi Systems BSBE \ Immunoturbidimetry (ITA)	mg/dL	265	2.4%	9	223		
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	mg/dL	228	16.4%	5	<100		
Roche cobas c 501/502/701/702 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	239	2.4%	5	141		
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	mg/dL	203	4.7%	43	1255		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	226	4.6%	34	675		
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	mg/dL	240	5.9%	28	1126		
Biostec A15/A25 Biostec \ Immunonephelometry (INA)	mg/dL	214	5.1%	27	718		
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	273	5.8%	25	408		
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	217	5.2%	20	458		
Siemens ADVIA 1800/2400 Siemens \ Immunoturbidimetry (ITA)	mg/dL	231	1.2%	<5	<100		
Mindray BS series MAKER \ Immunoturbidimetry (ITA)	mg/dL	118	24.6%	<5	<100		
Siemens ADVIA Centaur CP Dedicated \ Immunoturbidimetry (ITA)	mg/dL	260	1.5%	<5	<100		
Abbott ARCHITECT c series BSBE \ OTHER	mg/dL	231	1.8%	<5	<100		
Beckman Immage 800 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	229	4.9%	17	493		
Siemens BN II Siemens \ OTHER	mg/dL	277	8.5%	16	431		
TOSHIBA TBA MedicalSystem \ Immunoturbidimetry (ITA)	mg/dL	192	7.5%	16	254		
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	mg/dL	248	8.9%	14	398		
Beckman Immage 800 Beckman Coulter \ OTHER	mg/dL	209	5.9%	14	309		
Roche Hitachi/Hitachi Systems MAKER \ Immunoturbidimetry (ITA)	mg/dL	202	8.5%	10	303		
<b>(IgE)免疫球蛋白 IgE</b>							
Roche cobas c 501/502/701/702 MedicalSystem \ Immunoturbidimetry (ITA)	IU/mL	41.9	10.8%	5	148		
Siemens BN II Siemens \ OTHER	IU/mL	52.1	8.1%	16	437		
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	39.7	9.3%	<5	<100		
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	38.2	9.9%	<5	<100		
<b>(IgG)免疫球蛋白 IgG</b>							
Roche Hitachi/Hitachi Systems BSBE \ Immunoturbidimetry (ITA)	mg/dL	1408	5.6%	9	226		
Roche Hitachi/Hitachi Systems MAKER \ Immunoturbidimetry (ITA)	mg/dL	1459	16.2%	7	217		
Roche Hitachi/Hitachi Systems DongOu \ Immunoturbidimetry (ITA)	mg/dL	905	2.2%	7	132		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	1303	4.2%	50	1024		
Roche cobas c 501/502/701/702 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	1344	2.1%	5	141		
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	mg/dL	1213	5.6%	43	1249		
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	mg/dL	1234	5%	<5	<100		
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	mg/dL	1276	3.3%	28	1123		
Biostec A15/A25 Biostec \ Immunonephelometry (INA)	mg/dL	1252	6.5%	27	718		

<b>REF</b>	521-2	6x3mL	<b>EXP</b>	2024-04-04	<b>LOT</b>	31B51F8	更新日期:	2024-08
								4 / 6

		批号: 31B51F8					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数			
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	1410	8.5%	24	301		
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	1203	6.7%	20	460		
Siemens ADVIA 1800/2400 Siemens \ Immunoturbidimetry (ITA)	mg/dL	1235	1.3%	<5	<100		
Mindray BS series MAKER \ Immunoturbidimetry (ITA)	mg/dL	745	5.5%	<5	<100		
Siemens ADVIA Centaur CP Dedicated \ Immunoturbidimetry (ITA)	mg/dL	1325	3.6%	<5	<100		
Abbott ARCHITECT c series BSBE \ OTHER	mg/dL	1399	3.1%	<5	<100		
Beckman Immage 800 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	1311	4.5%	17	496		
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	mg/dL	1320	4.1%	16	453		
Siemens BN II Siemens \ OTHER	mg/dL	1409	8.2%	16	437		
TOSHIBA TBA MedicalSystem \ Immunoturbidimetry (ITA)	mg/dL	1208	14.2%	15	163		
<b>(IgM)免疫球蛋白 IgM</b>							
Roche Hitachi/Hitachi Systems BSBE \ Immunoturbidimetry (ITA)	mg/dL	176	4.2%	9	225		
Roche cobas 4000/6000/8000 Roche \ Immunonephelometry (INA)	mg/dL	162	9.4%	9	213		
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	mg/dL	203	3.1%	8	223		
Roche cobas c 501/502/701/702 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	164	2.8%	5	141		
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	mg/dL	165	3.4%	43	1252		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	165	3%	34	679		
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	mg/dL	168	6.6%	<5	<100		
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	mg/dL	160	2.3%	28	1117		
Biostec A15/A25 Biostec \ Immunonephelometry (INA)	mg/dL	128	6.9%	27	718		
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	156	6.9%	23	383		
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	162	6.6%	20	460		
Siemens ADVIA 1800/2400 Siemens \ Immunoturbidimetry (ITA)	mg/dL	170	1.6%	<5	<100		
Siemens ADVIA Centaur CP Dedicated \ Immunoturbidimetry (ITA)	mg/dL	172	1.7%	<5	<100		
Mindray BS series MAKER \ Immunoturbidimetry (ITA)	mg/dL	78.2	0.7%	<5	<100		
Abbott ARCHITECT c series BSBE \ OTHER	mg/dL	196	1.3%	<5	<100		
Beckman Immage 800 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	170	2.7%	17	497		
Siemens BN II Siemens \ OTHER	mg/dL	184	8%	16	438		
TOSHIBA TBA MedicalSystem \ Immunoturbidimetry (ITA)	mg/dL	143	5.4%	16	218		
Beckman Immage 800 Beckman Coulter \ OTHER	mg/dL	165	4.7%	14	306		
<b>(Kappa Light Chain)κ轻链</b>							
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	g/L	9.39	5.7%	43	1255		
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	g/L	9.56	7.2%	<5	<100		
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	g/L	9.34	10.3%	<5	<100		
Siemens BN II Dedicated \ Chemiluminescence Immunoassay (CLIA)	g/L	9.34	10.5%	<5	<100		
<b>(Lambda Light Chain)λ轻链</b>							
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	g/L	4.61	5.4%	43	1255		
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	g/L	4.57	4.6%	<5	<100		
Siemens BN II Dedicated \ Chemiluminescence Immunoassay (CLIA)	g/L	4.32	11.1%	<5	<100		
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	g/L	4.78	10.7%	<5	<100		
<b>(MA)微量白蛋白</b>							
Roche cobas c 501/502/701/702 Roche \ Immunoturbidimetry (ITA)	mg/L	34.9	8%	5	142		
<b>(PAB)前白蛋白</b>							
Roche Hitachi/Hitachi Systems CONDOR-TECO \ Immunoturbidimetry (ITA)	g/L	0.333	9.6%	6	174		
Roche cobas 4000/6000/8000 Roche \ Immunonephelometry (INA)	g/L	0.317	2.8%	5	<100		
Roche Hitachi/Hitachi Systems CONDOR-TECO \ Immunoturbidimetry (ITA)	%	36.6	6.6%	5	134		
Roche Hitachi/Hitachi Systems BSBE \ Immunoturbidimetry (ITA)	g/L	0.327	1.8%	34	629		

<b>REF</b>	521-2	6x3mL	<b>EXP</b>	2024-04-04	<b>LOT</b>	31B51F8	更新日期:	2024-08
								5 / 6

		批号: 31B51F8					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数			
Roche Hitachi/Hitachi Systems MedicalSystem \ Immunoturbidimetry (ITA)	g/L	0.315	5%	<5	<100		
Siemens ADVIA 1800/2400 Szybio \ Immunoturbidimetry (ITA)	g/L	0.315	2.9%	<5	115		
Beckman AU series HuaChenBio \ Immunonephelometry (INA)	g/L	0.323	3.7%	14	<100		
Beckman AU series KINGSBIO \ Immunoturbidimetry (ITA)	g/L	0.318	2.2%	13	369		
Siemens ADVIA 1800/2400 Siemens \ Immunoturbidimetry (ITA)	g/L	0.327	4.5%	<5	<100		
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA )	g/L	0.377	9%	<5	<100		
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	g/L	0.265	10.6%	<5	<100		
NIHON KOHDEN MEK Series MAKER \ Immunoturbidimetry (ITA)	g/L	0.294	2.6%	<5	<100		
<b>(RBP)视黄醇结合蛋白</b>							
Siemens ADVIA 1800/2400 BIOSINO \ Immunoturbidimetry (ITA)	mg/dL	4.9	11.8%	30	1235		
Roche Hitachi/Hitachi Systems MedicalSystem \ Immunoturbidimetry (ITA)	mg/dL	4.4	4.2%	<5	<100		
Siemens BN II Dedicated \ Chemiluminescence Immunoassay (CLIA )	mg/dL	4.43	10.4%	<5	<100		
Beckman AU series chongqiong zzybio \ Immunoturbidimetry (ITA)	mg/dL	4.73	2.3%	<5	<100		
<b>(RF)类风湿因子</b>							
Roche Hitachi/Hitachi Systems BSBE \ Immunoturbidimetry (ITA)	IU/mL	85	1.9%	9	226		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	IU/mL	59.6	7.9%	50	1032		
Beckman AU series HuaChenBio \ Immunonephelometry (INA)	IU/mL	67.3	1.5%	5	<100		
Roche cobas c 501/502/701/702 Dedicated \ Immunoturbidimetry (ITA)	IU/mL	65.7	1.8%	5	141		
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	IU/mL	57.8	7.3%	41	1224		
Roche cobas c 501/502/701/702 Roche \ Immunoturbidimetry (ITA)	IU/mL	81.1	2.5%	<5	<100		
Abbott ARCHITECT c series BSBE \ OTHER	IU/mL	97.1	5.2%	<5	<100		
Roche Hitachi/Hitachi Systems Wako \ Immunoturbidimetry (ITA)	IU/mL	58.8	9.1%	<5	<100		
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	IU/mL	60.3	4.5%	28	1117		
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	IU/mL	55.6	15.3%	25	310		
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	IU/mL	60	4%	21	549		
Roche cobas c 501/502/701/702 Roche \ Immunonephelometry (INA)	IU/mL	32.5	24.8%	<5	<100		
Mindray BS series Whitman \ Immunoturbidimetry (ITA)	IU/mL	38.2	1.1%	<5	<100		
Siemens ADVIA Centaur CP Dedicated \ Immunoturbidimetry (ITA)	IU/mL	74.5	1.8%	<5	<100		
Siemens ADVIA 1800/2400 Siemens \ Immunoturbidimetry (ITA)	IU/mL	71.1	1.5%	<5	134		
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	IU/mL	41.9	4.1%	16	459		
Siemens BN II Siemens \ OTHER	IU/mL	56	8.7%	16	423		
Beckman Immage 800 Beckman Coulter \ OTHER	IU/mL	58.5	7.4%	14	296		
TOSHIBA TBA BSBE \ Immunoturbidimetry (ITA)	IU/mL	90.9	5.4%	14	293		
<b>(TRF)转铁蛋白</b>							
Roche cobas c 501/502/701/702 MedicalSystem \ Immunoturbidimetry (ITA)	mg/dL	171	9.4%	5	142		
Abbott ARCHITECT c series BSBE \ OTHER	mg/dL	289	1%	<5	<100		
Roche Hitachi/Hitachi Systems LEADMAN \ Immunoturbidimetry (ITA)	mg/dL	267	2.8%	17	338		
Siemens BN II Siemens \ OTHER	mg/dL	263	10.6%	16	439		
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	mg/dL	246	5.4%	14	554		
Mindray BC-2000 series Mindray \ Immunoturbidimetry (ITA)	mg/dL	254	1%	<5	<100		
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA )	mg/dL	235	9.8%	<5	<100		
<b>(α-1-AG)α-1-酸性糖蛋白</b>							
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	g/L	0.865	9.9%	<5	<100		
<b>(α-1-AT)α-1-抗胰蛋白酶</b>							
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	123	11.4%	<5	<100		
<b>(α-1-G)α-1-球蛋白</b>							
TOSHIBA TBA LEADMAN \ Capillary Method	%	22.2	2.2%	<5	103		

<b>REF</b>	521-2	6x3mL	<b>EXP</b>	2024-04-04	<b>LOT</b>	31B51F8	更新日期：	2024-08 6 / 6
------------	-------	-------	------------	------------	------------	---------	-------	------------------

批号： 31B51F8					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
<b>(α1M)α1微球蛋白</b>					
Roche cobas c 501/502/701/702 Roche \ immunoturbidimetric	mg/L	22.3	15.4%	<5	<100
<b>(β-2-G)β-2-球蛋白</b>					
TOSHIBA TBA LEADMAN \ Capillary Method	%	3.94	4.3%	6	110
<b>(β2M)β-2-微球蛋白</b>					
Roche Hitachi/Hitachi Systems MedicalSystem \ Immunoturbidimetry (ITA)	ng/mL	3390	3.3%	9	239
Roche Hitachi/Hitachi Systems BSBE \ Immunoturbidimetry (ITA)	ng/mL	3390	2.7%	34	629
Abbott ARCHITECT c series BSBE \ OTHER	ng/mL	3654	5%	22	454
Mindray BS series Mindray \ OTHER	ng/mL	3904	5.5%	20	332
TOSHIBA TBA MedicalSystem \ OTHER	ng/mL	3470	9.1%	19	256
Beckman AU series chongqiong zybio \ Immunoturbidimetry (ITA)	ng/mL	3556	3%	14	391
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	ng/mL	3966	4.5%	12	234
NIHON KOHDEN MEK Series MAKER \ Immunoturbidimetry (ITA)	ng/mL	3809	1.9%	<5	<100
Beckman AU series HuaChenBio \ Immunonephelometry (INA)	ng/mL	4025	1.1%	<5	<100
Beckman Image 800 Dedicated \ Chemiluminescence Immunoassay (CLIA )	ng/mL	4300	14%	<5	<100
Siemens BN II Dedicated \ Chemiluminescence Immunoassay (CLIA )	ng/mL	3600	11.1%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Immunonephelometry (INA)	ng/mL	3655	0%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。