

<b>REF</b>	521-2	6x3mL	<b>EXP</b>	2026-01-26	<b>LOT</b>	31B52BE	更新日期:	2024-08
								1 / 3

		批号: 31B52BE					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数			
<b>(Alb)白蛋白</b>							
Roche Hitachi/Hitachi Systems Hitachi \ Immunoturbidimetry (ITA)	g/L	36.3	9.8%	<5	<100		
<b>(ASO)抗链球菌溶血素 O 抗体</b>							
Roche Hitachi/Hitachi Systems Whitman \ Immunoturbidimetry (ITA)	IU/mL	110	1.9%	9	470		
Mindray BS series Whitman \ Immunoturbidimetry (ITA)	IU/mL	118	3.7%	9	400		
Mindray BS series MedicalSystem \ OTHER	IU/mL	297	2%	8	382		
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	199	6.5%	<5	214		
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	231	2.9%	<5	102		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	IU/mL	172	11.5%	26	1410		
Siemens BN II Siemens \ Immunoturbidimetry (ITA)	IU/mL	258	10.3%	<5	<100		
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	IU/mL	186	1.6%	<5	<100		
<b>(C3)补体 C3</b>							
Roche Hitachi/Hitachi Systems MAKER \ Immunoturbidimetry (ITA)	mg/dL	93.3	5.8%	9	472		
Mindray BS series MedicalSystem \ OTHER	mg/dL	138	2.6%	8	384		
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	140	8.4%	<5	198		
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	168	4.7%	<5	102		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	145	3.8%	26	1410		
Mindray BS series MAKER \ Immunoturbidimetry (ITA)	mg/dL	84.4	6.6%	10	434		
Siemens BN II Siemens \ Immunoturbidimetry (ITA)	mg/dL	153	12.3%	<5	<100		
Beckman Immage 800 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	150	2.4%	<5	<100		
<b>(C4)补体 C4</b>							
Roche Hitachi/Hitachi Systems MAKER \ Immunoturbidimetry (ITA)	mg/dL	19.5	6.1%	9	472		
Mindray BS series MedicalSystem \ OTHER	mg/dL	34.7	1.6%	8	386		
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	38.3	7.8%	<5	212		
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	42.4	5.3%	<5	102		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	38.5	4.1%	26	1412		
Mindray BS series MAKER \ Immunoturbidimetry (ITA)	mg/dL	19.7	6.8%	10	458		
Siemens BN II Siemens \ Immunoturbidimetry (ITA)	mg/dL	48	10.1%	<5	<100		
Beckman Immage 800 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	40.5	3.5%	<5	<100		
<b>(Cer)铜蓝蛋白</b>							
Siemens ADVIA 1800/2400 BIOSINO \ Immunoturbidimetry (ITA)	g/L	32.4	1.7%	<5	<100		
<b>(CRP)C反应蛋白</b>							
Mindray BS series MedicalSystem \ OTHER	mg/dL	2.39	4.1%	8	191		
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	2.49	3.3%	<5	<100		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	2.44	2.9%	26	706		
Siemens BN II Siemens \ Immunoturbidimetry (ITA)	mg/dL	2.56	10.9%	<5	<100		
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	mg/dL	2.45	0.9%	<5	<100		
<b>(Hap (HPT) )触珠蛋白 (结合珠蛋白)</b>							
Roche Hitachi/Hitachi Systems Hitachi \ Immunoturbidimetry (ITA)	mg/L	1240	12.4%	<5	<100		
<b>(IgA )免疫球蛋白 IgA</b>							
Roche Hitachi/Hitachi Systems MAKER \ Immunoturbidimetry (ITA)	mg/dL	180	22.2%	9	466		
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	222	9.7%	<5	192		
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	315	4.3%	<5	102		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	245	6.1%	26	1402		
Mindray BS series MAKER \ Immunoturbidimetry (ITA)	mg/dL	112	11.6%	10	428		

<b>REF</b>	521-2	6x3mL	<b>EXP</b>	2026-01-26	<b>LOT</b>	31B52BE	更新日期:	2024-08
								2 / 3

		批号: 31B52BE					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数			
Siemens BN II Siemens \ Immunoturbidimetry (ITA)	mg/dL	265	12.4%	<5	<100		
Beckman Immage 800 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	259	2.1%	<5	<100		
<b>(IgG)免疫球蛋白 IgG</b>							
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	1390	6.1%	<5	210		
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	1517	5.9%	<5	102		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	1396	3.9%	26	1414		
Siemens BN II Siemens \ Immunoturbidimetry (ITA)	mg/dL	1520	12.3%	<5	<100		
Beckman Immage 800 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	1502	3.9%	<5	<100		
<b>(IgM)免疫球蛋白 IgM</b>							
Mindray BS series MAKER \ Immunoturbidimetry (ITA)	mg/dL	85.2	5.2%	9	404		
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	149	6.1%	<5	214		
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	175	3.6%	<5	102		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	150	3%	26	1412		
Roche Hitachi/Hitachi Systems MAKER \ Immunoturbidimetry (ITA)	mg/dL	92.4	11.8%	10	516		
Siemens BN II Siemens \ Immunoturbidimetry (ITA)	mg/dL	177	11.5%	<5	<100		
Beckman Immage 800 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	154	3.8%	<5	<100		
<b>(Kappa Light Chain)κ轻链</b>							
Siemens BN II Dedicated \ Chemiluminescence Immunoassay (CLIA )	g/L	4.01	5.7%	<5	<100		
Siemens BN II Siemens \ Immunoturbidimetry (ITA)	g/L	3.77	11.9%	<5	<100		
<b>(Lambda Light Chain)λ轻链</b>							
Siemens BN II Dedicated \ Chemiluminescence Immunoassay (CLIA )	g/L	2.3	6.3%	<5	<100		
Siemens BN II Siemens \ Immunoturbidimetry (ITA)	g/L	2.38	11.8%	<5	<100		
<b>(PAB)前白蛋白</b>							
Mindray BS series OTHER \ OTHER	g/L	0.387	3.1%	8	380		
Roche Hitachi/Hitachi Systems BSBE \ Immunoturbidimetry (ITA)	g/L	0.351	2.9%	26	1344		
Roche Hitachi/Hitachi Systems Hitachi \ Immunoturbidimetry (ITA)	g/L	0.36	12.9%	<5	<100		
<b>(RBP)视黄醇结合蛋白</b>							
Roche Hitachi/Hitachi Systems Hitachi \ Immunoturbidimetry (ITA)	mg/dL	6.27	11.9%	<5	<100		
<b>(RF)类风湿因子</b>							
Roche Hitachi/Hitachi Systems Whitman \ Immunoturbidimetry (ITA)	IU/mL	43	4.2%	9	472		
Mindray BS series MedicalSystem \ OTHER	IU/mL	85.4	2.8%	8	380		
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	59.6	5.1%	<5	102		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	IU/mL	80.2	6.5%	26	1422		
Mindray BS series Whitman \ Immunoturbidimetry (ITA)	IU/mL	43.8	5.9%	10	434		
Siemens BN II Siemens \ Immunoturbidimetry (ITA)	IU/mL	66.4	11.9%	<5	<100		
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	IU/mL	72.4	3.1%	<5	<100		
<b>(sTfR)可溶性转铁蛋白受体</b>							
Siemens BN II Siemens \ Immunoturbidimetry (ITA)	mg/L	0.96	11.4%	<5	<100		
<b>(TRF)转铁蛋白</b>							
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	257	6.4%	<5	<100		
Roche Hitachi/Hitachi Systems LEADMAN \ Immunoturbidimetry (ITA)	mg/dL	276	4.1%	14	387		
Roche Hitachi/Hitachi Systems Hitachi \ Immunoturbidimetry (ITA)	mg/dL	247	9.7%	<5	<100		
<b>(α-1-AG)α-1-酸性糖蛋白</b>							
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	g/L	1.01	8.9%	<5	<100		
Roche Hitachi/Hitachi Systems Hitachi \ Immunoturbidimetry (ITA)	g/L	0.68	9.6%	<5	<100		

<b>REF</b>	521-2	6x3mL	<b>EXP</b>	2026-01-26	<b>LOT</b>	31B52BE	更新日期：	2024-08 3 / 3
------------	-------	-------	------------	------------	------------	---------	-------	------------------

		批号：31B52BE			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
<b>(α-1-AT)α-1-抗胰蛋白酶</b>					
Roche Hitachi/Hitachi Systems Hitachi \ Immunoturbidimetry (ITA)	mg/dL	148	9%	<5	<100
<b>(β2M)β-2-微球蛋白</b>					
Mindray BS series MedicalSystem \ OTHER	ng/mL	4003	3.1%	8	380
Roche Hitachi/Hitachi Systems MedicalSystem \ Immunoturbidimetry (ITA)	ng/mL	3786	1.9%	<5	212
Roche Hitachi/Hitachi Systems BSBE \ Immunoturbidimetry (ITA)	ng/mL	3579	2.1%	26	1346
Siemens BN II Siemens \ Immunoturbidimetry (ITA)	ng/mL	3480	12%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。