

REF	521-2	6x3mL	EXP	2026-05-22	LOT	31B52C1	更新日期:	2024-08
								1 / 3

		批号: 31B52C1					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数			
(Alb)白蛋白							
Olympus AU400							
Olympus \ Bromocresol green	g/L	0.038	10.1%	<5	<100		
(ASO)抗链球菌溶血素 O 抗体							
Mindray BS series	IU/mL	277	5.5%	6	135		
MedicalSystem \ OTHER							
Abbott ARCHITECT c series	IU/mL	219	1.7%	5	102		
BSBE \ OTHER							
Siemens BN II	IU/mL	213	2.1%	<5	<100		
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)							
Beckman AU series	IU/mL	238	22%	<5	<100		
BSBE \ OTHER							
Mindray BS series	IU/mL	272	10%	<5	<100		
Mindray \ Immunoturbidimetry (ITA)							
Siemens BN ProSpec	IU/mL	242	1.3%	<5	<100		
Siemens \ Immunonephelometry (INA)							
Beckman Immage 800	IU/mL	166	4.8%	<5	<100		
Beckman Coulter \ Immunonephelometry (INA)							
KHB Polaris i2400	IU/mL	335	10%	<5	<100		
KHB \ Chemiluminescence Immunoassay (CLIA)							
Olympus AU400	IU/mL	265	10%	<5	<100		
Olympus \ Immunonephelometry (INA)							
Siemens ADVIA Centaur XP/XPT	IU/mL	251	0.4%	<5	<100		
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)							
Beckman AU series	IU/mL	264	10%	<5	<100		
Beckman Coulter \ Immunonephelometry (INA)							
(AT-III)抗凝血酶III							
Beckman Immage 800	g/L	0.236	10%	<5	<100		
Beckman Coulter \ Nephelometry							
(C3)补体 C3							
Mindray BS series	mg/dL	129	2.7%	6	135		
MedicalSystem \ OTHER							
Siemens BN II	mg/dL	152	4.3%	<5	<100		
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)							
KHB Polaris i2400	mg/dL	144	10%	<5	<100		
KHB \ Chemiluminescence Immunoassay (CLIA)							
Mindray BS series	mg/dL	144	10%	<5	<100		
Mindray \ Immunoturbidimetry (ITA)							
Beckman AU series	mg/dL	142	10%	<5	<100		
Beckman Coulter \ Immunonephelometry (INA)							
Beckman Immage 800	mg/dL	136	3.5%	<5	<100		
Beckman Coulter \ Immunonephelometry (INA)							
Siemens ADVIA Centaur XP/XPT	mg/dL	1.47	4.2%	<5	<100		
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)							
Olympus AU400	mg/dL	134	10%	<5	<100		
Olympus \ Immunonephelometry (INA)							
(C4)补体 C4							
Mindray BS series	mg/dL	32.5	1.3%	6	134		
MedicalSystem \ OTHER							
Siemens BN II	mg/dL	39.3	4.7%	<5	<100		
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)							
Beckman Immage 800	mg/dL	39.8	4.3%	<5	<100		
Beckman Coulter \ Immunonephelometry (INA)							
Beckman AU series	mg/dL	39	10%	<5	<100		
Beckman Coulter \ Immunonephelometry (INA)							
Siemens ADVIA Centaur XP/XPT	mg/dL	0.342	3.9%	<5	<100		
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)							
Mindray BS series	mg/dL	33	10%	<5	<100		
Mindray \ Immunoturbidimetry (ITA)							
Olympus AU400	mg/dL	39.3	10%	<5	<100		
Olympus \ Immunonephelometry (INA)							
KHB Polaris i2400	mg/dL	34	10%	<5	<100		
KHB \ Chemiluminescence Immunoassay (CLIA)							
(Cer)铜蓝蛋白							
Olympus AU400	g/L	0.483	10%	<5	<100		
Olympus \ Immunonephelometry (INA)							
(CRP)C反应蛋白							
Mindray BS series	mg/dL	2.82	17.3%	6	135		
MedicalSystem \ OTHER							
Siemens BN II	mg/dL	2.8	3%	<5	<100		
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)							
Mindray BS series	mg/dL	3.01	10%	<5	<100		
Mindray \ Immunoturbidimetry (ITA)							
Olympus AU400	mg/dL	2.98	9.7%	<5	<100		
Olympus \ Immunonephelometry (INA)							
TOSHIBA TBA	mg/dL	2.2	6.1%	<5	<100		
TOSHIBA \ Immunoturbidimetry (ITA)							
Siemens BN ProSpec	mg/dL	2.76	0.4%	<5	<100		
Siemens \ Immunonephelometry (INA)							
Beckman AU series	mg/dL	2.88	1.8%	<5	<100		
Beckman Coulter \ Immunonephelometry (INA)							

REF	521-2	6x3mL	EXP	2026-05-22	LOT	31B52C1	更新日期:	2024-08
								2 / 3

批号: 31B52C1						
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数		
(Ferritin)铁蛋白						
Siemens BN II						
Siemens \ Immunoturbidimetry (ITA)	ng/mL	116	10%	<5	<100	
Olympus AU400						
Olympus \ Immunoturbidimetry (ITA)	ng/mL	135	10%	<5	<100	
(Hap (HPT))触珠蛋白 (结合珠蛋白)						
Abbott ARCHITECT c series						
chongqiong zybio \ OTHER	mg/L	1321	5%	9	223	
Beckman AU series						
chongqiong zybio \ OTHER	mg/L	1422	10.7%	<5	<100	
Olympus AU400						
Olympus \ Immunonephelometry (INA)	mg/L	1350	10%	<5	<100	
(IgA)免疫球蛋白 IgA						
Siemens BN II						
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	326	2.4%	<5	<100	
Siemens ADVIA Centaur XP/XPT						
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	mg/dL	262	2.5%	<5	<100	
Mindray BS series						
Mindray \ Immunoturbidimetry (ITA)	mg/dL	293	10%	<5	<100	
Beckman AU series						
Beckman Coulter \ Immunonephelometry (INA)	mg/dL	292	10%	<5	<100	
Olympus AU400						
Olympus \ Immunonephelometry (INA)	mg/dL	271	10%	<5	<100	
Abbott ARCHITECT c series						
DiaSys \ Immunoturbidimetry (ITA)	mg/dL	267	10%	<5	<100	
KHB Polaris i2400						
KHB \ Chemiluminescence Immunoassay (CLIA)	mg/dL	326	10%	<5	<100	
(IgE)免疫球蛋白 IgE						
Siemens BN II						
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	48.5	2.6%	<5	<100	
Siemens BN ProSpec						
Siemens \ Immunonephelometry (INA)	IU/mL	48.8	4%	<5	<100	
Beckman Immage 800						
Beckman Coulter \ Immunonephelometry (INA)	IU/mL	47.1	10%	<5	<100	
(IgG)免疫球蛋白 IgG						
Siemens BN II						
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	1564	4.5%	<5	<100	
Siemens ADVIA Centaur XP/XPT						
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	mg/dL	1607	2%	<5	<100	
Beckman AU series						
Beckman Coulter \ Immunonephelometry (INA)	mg/dL	1471	10%	<5	<100	
Abbott ARCHITECT c series						
DiaSys \ Immunoturbidimetry (ITA)	mg/dL	1415	10%	<5	<100	
Mindray BS series						
Mindray \ Immunoturbidimetry (ITA)	mg/dL	1240	10%	<5	<100	
Olympus AU400						
Olympus \ Immunonephelometry (INA)	mg/dL	1380	10%	<5	<100	
KHB Polaris i2400						
KHB \ Chemiluminescence Immunoassay (CLIA)	mg/dL	1520	10%	<5	<100	
Beckman Immage 800						
Beckman Coulter \ Immunonephelometry (INA)	mg/dL	1411	6.3%	<5	<100	
(IgM)免疫球蛋白 IgM						
Siemens BN II						
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	210	2.9%	<5	<100	
Beckman Immage 800						
Beckman Coulter \ Immunonephelometry (INA)	mg/dL	182	4.2%	<5	<100	
Mindray BS series						
Mindray \ Immunoturbidimetry (ITA)	mg/dL	172	10%	<5	<100	
Siemens ADVIA Centaur XP/XPT						
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	mg/dL	161	1.6%	<5	<100	
Abbott ARCHITECT c series						
DiaSys \ Immunoturbidimetry (ITA)	mg/dL	163	10%	<5	<100	
Olympus AU400						
Olympus \ Immunonephelometry (INA)	mg/dL	165	10%	<5	<100	
Beckman AU series						
Beckman Coulter \ Immunonephelometry (INA)	mg/dL	179	10%	<5	<100	
KHB Polaris i2400						
KHB \ Chemiluminescence Immunoassay (CLIA)	mg/dL	262	10%	<5	<100	
(Kappa Light Chain)κ轻链						
Beckman Immage 800						
Beckman Coulter \ Immunonephelometry (INA)	g/L	10.6	10%	<5	<100	
(Lambda Light Chain)λ 轻链						
Beckman Immage 800						
Beckman Coulter \ Immunonephelometry (INA)	g/L	5.68	10%	<5	<100	
(PAB)前白蛋白						
Mindray BS series						
OTHER \ OTHER	g/L	0.318	2.6%	6	135	
Beckman AU series						
BSBE \ OTHER	g/L	3.46	18.4%	<5	<100	
Abbott ARCHITECT c series						
BSBE \ OTHER	g/L	3.09	4.9%	10	218	

REF	521-2	6x3mL	EXP	2026-05-22	LOT	31B52C1	更新日期:	2024-08
								3 / 3

		批号: 31B52C1			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
TOSHIBA TBA TOSHIBA \ Immunoturbidimetry (ITA)	g/L	0.357	1.1%	<5	<100
KHB Polaris i2400 KHB \ Chemiluminescence Immunoassay (CLIA)	g/L	0.321	10%	<5	<100
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	g/L	0.32	10%	<5	<100
Olympus AU400 Olympus \ Immunonephelometry (INA)	g/L	0.315	10%	<5	<100
Siemens ADVIA Centaur XP/XPT Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	g/L	287	0.4%	<5	<100
(RBP)视黄醇结合蛋白					
Siemens BN II Siemens \ Immunonephelometry (INA)	mg/dL	4.73	9.9%	<5	<100
(RF)类风湿因子					
Mindray BS series MedicalSystem \ OTHER	IU/mL	78.4	2.4%	6	118
Abbott ARCHITECT c series BSBE \ OTHER	IU/mL	77.2	2.7%	5	105
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	IU/mL	71.6	10.6%	<5	<100
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	54.1	6.5%	<5	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	IU/mL	67.1	10%	<5	<100
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	IU/mL	55.2	3.3%	<5	<100
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	IU/mL	60.1	10%	<5	<100
TOSHIBA TBA TOSHIBA \ Immunoturbidimetry (ITA)	IU/mL	74.2	0.7%	<5	<100
KHB Polaris i2400 KHB \ Chemiluminescence Immunoassay (CLIA)	IU/mL	76.9	10%	<5	<100
(TRF)转铁蛋白					
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	237	2.8%	<5	<100
Olympus AU400 Olympus \ Immunoturbidimetry (ITA)	mg/dL	235	10%	<5	<100
(α-1-AG)α-1-酸性糖蛋白					
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	g/L	1.1	7.9%	<5	<100
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	g/L	0.698	1.7%	<5	<100
Olympus AU400 Olympus \ Immunoturbidimetry (ITA)	g/L	0.772	10%	<5	<100
(α-1-AT)α-1-抗胰蛋白酶					
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	152	5.4%	<5	<100
Olympus AU400 Olympus \ Immunoturbidimetry (ITA)	mg/dL	149	10%	<5	<100
(α2M)α-2-巨球蛋白					
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	mg/L	2180	10%	<5	<100
(β2M)β-2-微球蛋白					
Mindray BS series MedicalSystem \ OTHER	ng/mL	4189	1.4%	6	135
Beckman AU series BSBE \ OTHER	ng/mL	3596	14.2%	<5	<100
Abbott ARCHITECT c series BSBE \ OTHER	ng/mL	3621	5.6%	10	217
Olympus AU400 Olympus \ Immunonephelometry (INA)	ng/mL	4280	9.3%	<5	<100
KHB Polaris i2400 KHB \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3340	10%	<5	<100
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	ng/mL	3970	10%	<5	<100
TOSHIBA TBA TOSHIBA \ Immunoturbidimetry (ITA)	ng/mL	2828	1%	<5	<100
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	ng/mL	3980	2.8%	<5	<100

注: 此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来, 仅供学习、参考之用。因所用技术、仪器和试剂的不同, 或因制造商检测方法的改变, 均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求, 实验室须遵循相关技术规范确立自己的均值和可接受范围。