

REF	521-3	6x3mL	EXP	2024-04-03	LOT	31B7908	更新日期:	2024-08
								1 / 3

		批号: 31B7908					
项目\仪器\试剂方法	单位	均值\CV		Lab\点数			
(Alb)白蛋白							
Roche cobas 4000/6000/8000 Roche \ Bromocresol green	g/L	57.8	12%	<5	<100		
(ASO)抗链球菌溶血素 O 抗体							
Beckman Immage 800 Beckman Coulter \ OTHER	IU/mL	188	5.9%	8	176		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	IU/mL	205	6%	31	678		
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	IU/mL	206	4.3%	13	<100		
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	IU/mL	199	4.3%	10	292		
(AT-III)抗凝血酶III							
Stago STA series Dedicated \ Nephelometry	%	48	12%	<5	<100		
(C3)补体 C3							
Beckman Immage 800 Beckman Coulter \ OTHER	mg/dL	205	5.1%	8	191		
Roche Hitachi/Hitachi Systems sysmex \ Immunoturbidimetry (ITA)	mg/dL	227	10.7%	5	<100		
Roche cobas c 501/502/701/702 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	231	3.2%	5	143		
TOSHIBA TBA MedicalSystem \ Immunoturbidimetry (ITA)	mg/dL	227	3%	<5	<100		
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	247	2.9%	<5	<100		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	203	4.5%	31	662		
Abbott ARCHITECT c series BSBE \ OTHER	mg/dL	216	3%	<5	<100		
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	mg/dL	209	3.9%	17	<100		
Siemens BN II Siemens \ OTHER	mg/dL	244	7.9%	16	444		
Beckman Immage 800 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	213	3.6%	10	293		
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	mg/dL	211	1.3%	<5	<100		
(C4)补体 C4							
Beckman Immage 800 Beckman Coulter \ OTHER	mg/dL	55.4	5.3%	8	186		
Roche Hitachi/Hitachi Systems sysmex \ Immunoturbidimetry (ITA)	mg/dL	43.8	10.6%	7	<100		
Roche cobas c 501/502/701/702 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	48.4	4.9%	5	144		
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	56.4	3.5%	<5	<100		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	54.7	5.5%	31	675		
TOSHIBA TBA MedicalSystem \ Immunoturbidimetry (ITA)	mg/dL	48.9	2.6%	<5	<100		
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	mg/dL	44.9	2.3%	<5	<100		
Abbott ARCHITECT c series BSBE \ OTHER	mg/dL	48.5	1.6%	<5	<100		
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	mg/dL	55.4	4.5%	17	<100		
Siemens BN II Siemens \ OTHER	mg/dL	54.6	7.2%	16	424		
Beckman Immage 800 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	56.3	3.7%	10	291		
(Cer)铜蓝蛋白							
Siemens BN II Siemens \ OTHER	g/L	0.778	9.1%	16	429		
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	g/L	0.452	3.6%	10	<100		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	g/L	0.45	12%	<5	<100		
(CRP)C反应蛋白							
Beckman Immage 800 Beckman Coulter \ OTHER	mg/dL	7.53	3.6%	8	188		
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	7.49	5%	31	666		
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	mg/dL	7.71	3.5%	17	<100		
Siemens BN II Siemens \ OTHER	mg/dL	6.37	7%	16	427		
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	mg/dL	7.56	4.6%	10	294		
(Ferritin)铁蛋白							
Siemens BN II Siemens \ OTHER	ng/mL	289	8.6%	16	432		
Roche Hitachi/Hitachi Systems LEADMAN \ Immunoturbidimetry (ITA)	ng/mL	235	1.5%	15	318		
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	341	12%	<5	<100		

REF	521-3	6x3mL	EXP	2024-04-03	LOT	31B7908	更新日期:	2024-08
								2 / 3

		批号: 31B7908			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
(Hap (HPT))触珠蛋白 (结合珠蛋白)					
Siemens BN II	mg/L	2274	11.1%	16	442
Siemens \ OTHER					
(HbA1c (%))糖化血红蛋白A1c (%)					
Roche cobas c 501/502/701/702	%	10.2	1.5%	5	147
Dedicated \ immunoturbidimetric					
(hs-CRP)超敏C反应蛋白					
TOSHIBA TBA	mg/L	56.1	6%	<5	<100
OTHER \ OTHER					
Abbott ARCHITECT c series	mg/L	60	8.5%	14	272
BSBE \ OTHER					
(IgA)免疫球蛋白 IgA					
Roche Hitachi/Hitachi Systems	mg/dL	377	15%	8	<100
sysmex \ Immunoturbidimetry (ITA)					
Beckman Immage 800	mg/dL	338	4.3%	8	193
Beckman Coulter \ OTHER					
Roche cobas c 501/502/701/702	mg/dL	360	2.3%	5	142
Dedicated \ Immunoturbidimetry (ITA)					
Mindray BS series	mg/dL	388	3%	<5	<100
Mindray \ Immunoturbidimetry (ITA)					
Beckman Immage 800	mg/dL	341	4.5%	31	678
Beckman Coulter \ Immunonephelometry (INA)					
TOSHIBA TBA	mg/dL	334	2.1%	<5	<100
MedicalSystem \ Immunoturbidimetry (ITA)					
Abbott ARCHITECT c series	mg/dL	374	1.4%	<5	<100
BSBE \ OTHER					
Beckman Immage 800	mg/dL	330	3.6%	17	<100
OTHER \ Immunonephelometry (INA)					
Siemens BN II	mg/dL	447	7.4%	16	432
Siemens \ OTHER					
Beckman Immage 800	mg/dL	351	2.8%	10	291
Dedicated \ Immunoturbidimetry (ITA)					
TOSHIBA TBA	mg/dL	363	2%	<5	<100
OTHER \ Immunoturbidimetry (ITA)					
(IgE)免疫球蛋白 IgE					
Roche cobas c 501/502/701/702	IU/mL	59	8.5%	5	148
MedicalSystem \ Immunoturbidimetry (ITA)					
Siemens BN II	IU/mL	76.6	5.8%	16	437
Siemens \ OTHER					
Beckman Immage 800	IU/mL	59.8	12%	<5	<100
Beckman Coulter \ Immunonephelometry (INA)					
(IgG)免疫球蛋白 IgG					
Beckman Immage 800	mg/dL	1855	4.6%	9	189
Beckman Coulter \ OTHER					
Roche Hitachi/Hitachi Systems	mg/dL	1731	13.2%	6	<100
sysmex \ Immunoturbidimetry (ITA)					
Roche cobas c 501/502/701/702	mg/dL	2002	2.9%	5	141
Dedicated \ Immunoturbidimetry (ITA)					
Mindray BS series	mg/dL	2193	3.8%	<5	<100
Mindray \ Immunoturbidimetry (ITA)					
Beckman Immage 800	mg/dL	1887	4.4%	31	683
Beckman Coulter \ Immunonephelometry (INA)					
TOSHIBA TBA	mg/dL	2387	4%	<5	<100
MedicalSystem \ Immunoturbidimetry (ITA)					
Abbott ARCHITECT c series	mg/dL	2172	3.4%	<5	<100
BSBE \ OTHER					
Beckman Immage 800	mg/dL	1884	3.6%	17	<100
OTHER \ Immunonephelometry (INA)					
Siemens BN II	mg/dL	2102	7.1%	16	433
Siemens \ OTHER					
Beckman Immage 800	mg/dL	1926	3%	10	292
Dedicated \ Immunoturbidimetry (ITA)					
TOSHIBA TBA	mg/dL	1756	1.5%	<5	<100
OTHER \ Immunoturbidimetry (ITA)					
(IgM)免疫球蛋白 IgM					
Beckman Immage 800	mg/dL	255	4.8%	8	188
Beckman Coulter \ OTHER					
Roche Hitachi/Hitachi Systems	mg/dL	262	4.1%	5	<100
sysmex \ Immunoturbidimetry (ITA)					
Roche cobas c 501/502/701/702	mg/dL	253	2.7%	5	142
Dedicated \ Immunoturbidimetry (ITA)					
TOSHIBA TBA	mg/dL	232	1.3%	<5	<100
MedicalSystem \ Immunoturbidimetry (ITA)					
Beckman Immage 800	mg/dL	255	4.2%	31	677
Beckman Coulter \ Immunonephelometry (INA)					
Mindray BS series	mg/dL	175	1.9%	<5	<100
Mindray \ Immunoturbidimetry (ITA)					
Roche Hitachi/Hitachi Systems	mg/dL	265	5.5%	<5	<100
OTHER \ Immunoturbidimetry (ITA)					
Beckman Immage 800	mg/dL	255	2.6%	17	<100
OTHER \ Immunonephelometry (INA)					
Siemens BN II	mg/dL	293	6.2%	16	441
Siemens \ OTHER					
Beckman Immage 800	mg/dL	261	2.6%	10	293
Dedicated \ Immunoturbidimetry (ITA)					
TOSHIBA TBA	mg/dL	242	2.4%	<5	<100
OTHER \ Immunoturbidimetry (ITA)					

REF 521-3	6x3mL	EXP 2024-04-03	LOT 31B7908	更新日期: 2024-08
				3 / 3

批号: 31B7908					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
(Kappa Light Chain)κ轻链					
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	g/L	16.3	3.6%	17	<100
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	g/L	15.5	12%	<5	<100
(Lambda Light Chain)λ 轻链					
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	g/L	8.31	5.4%	17	<100
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	g/L	8.13	12%	<5	<100
(MA)微量白蛋白					
Roche cobas c 501/502/701/702 Roche \ Immunoturbidimetry (ITA)	mg/L	133	5%	5	142
(PAB)前白蛋白					
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	g/L	0.476	2%	7	<100
Roche Hitachi/Hitachi Systems BSBE \ Immunoturbidimetry (ITA)	g/L	0.491	1.8%	30	593
Roche cobas 4000/6000/8000 DiaSys \ Immunonephelometry (INA)	g/L	0.48	12%	<5	<100
(RBP)视黄醇结合蛋白					
Roche cobas 4000/6000/8000 WanTai \ Immunonephelometry (INA)	mg/dL	6.56	12%	<5	<100
(RF)类风湿因子					
Beckman Immage 800 Beckman Coulter \ OTHER	IU/mL	121	4%	8	142
TOSHIBA TBA OTHER \ OTHER	IU/mL	129	2.8%	6	<100
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	104	1.7%	6	<100
Roche cobas c 501/502/701/702 Dedicated \ Immunoturbidimetry (ITA)	IU/mL	97.2	2.7%	5	141
TOSHIBA TBA MedicalSystem \ Immunoturbidimetry (ITA)	IU/mL	136	4.2%	<5	<100
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	IU/mL	111	6.5%	31	667
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	IU/mL	89.1	4.8%	<5	<100
Abbott ARCHITECT c series BSBE \ OTHER	IU/mL	90.9	2.8%	<5	<100
Siemens BN II Siemens \ OTHER	IU/mL	86.4	7.9%	16	430
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	IU/mL	112	4.1%	13	<100
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	IU/mL	110	3.4%	10	293
(TRF)转铁蛋白					
Beckman Immage 800 Beckman Coulter \ OTHER	mg/dL	353	6.1%	9	152
Siemens BN II Siemens \ OTHER	mg/dL	398	9.9%	16	429
Roche Hitachi/Hitachi Systems LEADMAN \ Immunoturbidimetry (ITA)	mg/dL	381	2.9%	15	318
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	mg/dL	357	12%	<5	<100
(α-1-AT)α-1-抗胰蛋白酶					
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	201	10.4%	<5	<100
(α1M)α1微球蛋白					
Roche cobas c 501/502/701/702 Roche \ immunoturbidimetric	mg/L	90.7	4%	<5	<100
(β2M)β-2-微球蛋白					
TOSHIBA TBA OTHER \ OTHER	ng/mL	6332	1.8%	8	<100
TOSHIBA TBA MedicalSystem \ OTHER	ng/mL	4664	4.1%	<5	<100
Roche Hitachi/Hitachi Systems BSBE \ Immunoturbidimetry (ITA)	ng/mL	5017	2.6%	30	593
Roche Hitachi/Hitachi Systems OTHER \ OTHER	ng/mL	5422	4.8%	<5	<100
Mindray BS series Mindray \ OTHER	ng/mL	5962	5.5%	<5	<100
Abbott ARCHITECT c series BSBE \ OTHER	ng/mL	5174	4%	14	282
Roche cobas 4000/6000/8000 JS (Shanghai) \ Immunonephelometry (INA)	ng/mL	5125	12%	<5	<100

注: 本报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来, 仅供学习、参考之用。因所用技术、仪器和试剂的不同, 或因制造商检测方法的改变, 均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求, 实验室须遵循相关技术规范确立自己的均值和可接受范围。